

# E-21 SINGLE USE FLUID MANAGEMENT SYSTEM

A SINGLE USE ENDOSCOPIC PUMP can be used as a replacement for existing Arthroscopy pumps or as a backup.



## **Pump-tube bundle system**

no need for complicated tubing installation, just plug in and use

01



## **Battery operated (5 hours of battery operation)**

Hassle-free use without hanging power lines

02



## **Wide pressure range (20-100mmHg)**

applicable to different types of endoscopic surgery

03



## **Lower pressure detectable (from 20mmHg)**

suitable for spinal endoscopic surgeries that require a small amount of saline inflow

04